


Experiment(s):

1. Simulation of radioactive decay (exponential decay) using rolling of dice

(For more details, procedure & manual visit: www.kamaljeeth.net)

Reference : *Lab Experiments Journal vol-11, No.2, Page-146*



(a)

Experiment setup consists:

- a) Board
- b) Dice set

Specifications:

a) Board

Size: 500 mm x 500 mm
Material: Wood

b) Dice set

Quantity: 100 nos - 5 nos spare
Size: 15 mm x 15 mm x 15 mm



(b)



KAMALJEETH INSTRUMENTS

Address: No. 610, 5th main, 8th cross Tatanagar, Bangalore - 560092, INDIA

Website: www.kamaljeeth.net, Email: labexperiments@kamaljeeth.net

ESTD. 1990

3 years manufacturing
warranty